

Enhancing imaging contrast of absorptive object using a backscattering matrix


Hakseok Ko, Wonjun Choi, Moonseok Kim, Mooseok Jang, and Wonshik Choi

Center for Molecular Spectroscopy and Dynamics, Institute for Basic Science (IBS), Seoul 02841, Korea

Department of Physics, Korea University, Seoul 136-713 Korea

Introduction

Multiple scattering events restricts the range of imaging depth



Eigenchannels differentiates objects depending on their absorptivity


- Random matrix theory predicts perfect reflection
- Singular value decomposition of the backscattering matrix

$$t = U \tau V^\dagger$$

V : Unitary matrix connecting input waves into eigenchannels

U : Unitary matrix connecting eigenchannels into output channels

τ : Real non-negative diagonal matrix




O. N. Dorokhov, Solid State Comm. 51, No. 6, 381 (1984)
Choi, W. et al. Phys. Rev. B 83, 134207 (2011).


- Each eigenchannel takes proper spatial distribution depending on the absorptivity distribution of the scattering medium

Materials and method

Definition of backscattering matrix





Experimental setup



- 633 nm He-Ne laser
- Scattering medium: Mixture of PDMS and ZnO powder
- Absorptive target: 3 μm & 10 μm blue microsphere
- Experimental condition
 - Illumination area: 21.6 x 21.6 μm
 - Illumination NA: 0.6
 - Recording NA: 0.6
 - Recorded image: 1800

Results

Absorptive area on reflective mirror surface




- From the eigenvalue curve, reflective area (B) and absorptive area can be identified

- We could selectively image the area with specific absorptivity, by selectively summing the eigenchannel images identified from the eigenvalue curve



- Scale bar 3 μm

Absorptive area on scattering medium



- Mixture of PDMS and Zinc oxide was used as a scattering medium
- We also could selectively image the reflective area (scattering medium) and the absorptive area
- Open eigenchannel tends to illuminates the scattering medium to make high reflectivity
- Closed eigenchannel tends to illuminates the scattering medium to make low reflectivity
- Scale bar 3 μm

Comparison with incoherent illumination imaging



- We divided open eigenchannel image by closed eigenchannel (div)

- By this, we improved the imaging contrast of absorptive objects by ~20 % in comparison with the conventional incoherent illumination imaging

Conclusion

- We recorded the backscattering matrix of a scattering medium with absorptive targets
- We identified reflection eigenchannels of the scattering medium
- Eigenchannels identified from the singular value decomposition of the backscattering matrix differentiated objects depending on their absorptivity
- By using the property, we improved the imaging contrast of absorptive objects by ~20 % in comparison with the conventional incoherent illumination imaging

Reference

- [1] Moonseok Kim, Youngwoon Choi, Changhyeong Yoon, Wonjun Choi, Jaisoon Kim, Q-Han Park and Wonshik Choi, "Maximal energy transport through disordered media with the implementation of transmission eigenchannels," Nature Photonics, 6 581 (2012)
- [2] Moonseok Kim, Wonjun Choi, Changhyeong Yoon, Guang Hoon Kim, Seunghyun Kim, Gi-Ra Yi, Q-Han Park, and Wonshik Choi, "Exploring anti-reflection modes in disordered media," Optics express, 23, 12740-12749 (6 May 2015)
- [3] Wonjun Choi, Allard P. Mosk, Q-Han Park and Wonshik Choi, "Transmission eigenchannels in a disordered medium," Physical ReviewB 83, 134207 (2011)